

STORZ

KARL STORZ — ENDOSKOPE

en **Instructions for use**
C-MAC® Video Laryngoscopes Series 8403



04-2021

Copyright ©

All product illustrations, product descriptions, and texts are the intellectual property of KARL STORZ SE & Co. KG.

Their use and reproduction by third parties require the express approval of KARL STORZ SE & Co. KG.

All rights reserved.



Table of contents

1	General information.....	5
1.1	Read the instructions for use.....	5
1.2	Read the instructions for use of combinable products.....	5
1.3	Scope.....	5
1.4	Description of warning messages.....	6
2	Normal use.....	7
2.1	Intended use.....	7
2.2	Indications for use.....	7
2.3	Contraindications.....	7
2.4	Target user populations.....	7
2.5	Patient groups.....	7
3	Safety.....	8
3.1	Serious incidents.....	8
3.2	Correct handling.....	8
3.3	Unsterile product.....	8
3.4	Damaged products.....	8
3.5	Combination with other components.....	9
3.6	Electromagnetic interference.....	9
3.7	Patient leakage current.....	10
3.8	National guidelines on airway management.....	10
3.9	Hot components.....	10
3.10	High light intensity.....	10
4	Product description.....	11
4.1	Product overview.....	11
4.2	Possible combinations.....	12
4.3	Technical specifications.....	12
4.4	Symbols employed.....	13
4.4.1	Symbols on the packaging.....	13
4.4.2	Symbols on the product.....	14
4.5	Ambient conditions.....	14
5	Preparation.....	15
5.1	Unpacking the product.....	15
5.2	Testing the product.....	15
5.2.1	Visual inspection.....	15
5.2.2	Functional test.....	17
5.3	Connecting the C-MAC® Monitor.....	17
5.4	Connecting the C-MAC® Pocket Monitor.....	19
5.5	Connecting C-Hub® II.....	20
6	Application.....	21
6.1	Using the product.....	21
6.2	Inserting an oxygen or suction catheter.....	21
7	Disassembly.....	23
7.1	Disassembling the product.....	23
8	Maintenance, servicing, repairs, and disposal.....	24
8.1	Repairs to the product.....	24
8.2	Disposing of the product.....	24
9	Accessories and spare parts.....	25
9.1	Accessories.....	25
10	Electromagnetic compatibility.....	26
10.1	Electromagnetic compatibility.....	26

11 Fault correction 27
11.1 Product malfunctions 27
12 Subsidiaries..... 28

1 General information

1.1 Read the instructions for use

If the instructions for use are not followed, patients, users, and third parties may be injured or the product may be damaged.

- ▶ Read the instructions for use carefully and follow all the safety notes and warnings.
- ▶ Keep the instructions for use clearly visible next to the product.

1.2 Read the instructions for use of combinable products

If the instructions for use of combinable products are not followed, patients, users, and third parties may be injured or the product may be damaged.

- ▶ Read the instructions for use of the combinable products carefully and follow all the safety notes and warnings.

1.3 Scope

These instructions for use are valid for the following products:

Product	Item number
C-MAC® Video Laryngoscope	8403AX
	8403AXC
	8403BX
	8403BXC
	8403HX
	8403DXC
	8403EXC
	8403GXC
	8403KXC
	8403MXC
	8403HXP
	8403NXC
Video connecting cable	8403X

1.4 Description of warning messages

To prevent any injury to persons or damage to property, the warnings and safety notes in the instructions for use must be observed. The warning messages describe the following levels of danger.

▲ WARNING**WARNING**

Designates a possible imminent risk. If this is not avoided, it could lead to death or serious injuries.

▲ CAUTION**CAUTION**

Designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

NOTICE**ATTENTION**

Designates a possibly harmful situation. If this is not avoided, the products could be damaged.

2 Normal use

2.1 Intended use

Video laryngoscopes (CMOS) are used for visualizing the respiratory tracts and vocal cords during endotracheal intubation and for the inspection and examination of the upper respiratory tract. Video laryngoscopes (CMOS) are designed for transient use in invasive procedures through a body orifice.

C-MAC® Monitor accessories are used for transmitting energy. C-MAC® Monitor accessories are used in conjunction with the C-MAC® System for visualizing the respiratory tracts and vocal cords during endotracheal intubation and for inspecting and examining the upper respiratory tract. C-MAC® Monitor accessories are non-invasive and are designed for transient use in invasive procedures through a body orifice.

2.2 Indications for use

The use of laryngoscopes and video laryngoscopes is indicated if, in the opinion of the responsible physician, inspection of the upper respiratory tract or endotracheal intubation is indicated.

2.3 Contraindications

The use of laryngoscopes and video laryngoscopes is contraindicated if, in the opinion of the responsible physician, use is contraindicated or the patient is not able to undergo surgery or anesthesia due to his or her general condition. Laryngoscopes and video laryngoscopes must not be used for procedures in direct contact with the central nervous system (CNS) and central cardiovascular system.

2.4 Target user populations

The medical device may only be used by doctors and medical assistants with a relevant specialist qualification.

2.5 Patient groups

There are no restrictions in terms of patient groups for this product.

3 Safety

3.1 Serious incidents

According to the Medical Device Regulation (MDR), a “serious incident” includes incidents that directly or indirectly had, could have had, or could have any of the following consequences (MDR, Art. 2, No. 65 [1]):

- Death of a patient, user, or another person
- Temporary or permanent serious deterioration in the medical condition of a patient, user, or another person
- A serious threat to public health
- ▶ The manufacturer and appropriate authority must be notified of all serious incidents.

3.2 Correct handling

If the product is not handled correctly, patients, users, and third parties may be injured.

- ▶ Only persons with the necessary medical qualification and who are acquainted with the application of the product may work with it.
- ▶ Check that the product is suitable for the procedure prior to use.
- ▶ Check the product for the following points before and after every use:
 - Completeness
 - Good working order
 - Rough surfaces left inadvertently
 - Sharp corners
 - Burred edges
 - Correct assembly of the components
 - Functionality
- ▶ Do not leave broken-off components inside the patient.
- ▶ Do not overload the product with mechanical stress.
- ▶ Do not bend bent products back to their original position.

3.3 Unsterile product

The product is not sterile when delivered. The use of non-sterile products poses a risk of infection for patients, users, and third parties.

- ▶ Reprocess the product in line with the reprocessing instructions before initial use and every subsequent use.

3.4 Damaged products

Damaged products can result in injury to patients, users, or third parties.

- ▶ Before each use, check all components of the product for damage.
- ▶ Do not use damaged products.

3.5 Combination with other components

The use of unauthorized devices and components or unauthorized changes to the product can result in injuries.

Additional devices connected to electrical medical equipment must comply with the relevant IEC or ISO standards. Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3rd edition of IEC 60601-1).

- ▶ Only combine the product with devices and components that are approved for joint use by the manufacturer.
- ▶ Comply with national and local regulations.
- ▶ Observe the instruction manuals and interface specifications of the devices and components used in combination.
- ▶ Only use devices and components that have standardized interfaces and do not breach the intended use of the product.
- ▶ Only make changes to the product if these changes are approved by KARL STORZ.

3.6 Electromagnetic interference

Medical electrical products are subject to special precautions regarding electromagnetic compatibility and must be installed and commissioned according to the tables on electromagnetic compatibility. If other products (e.g. for MRT, CT, diathermy, electrocautery, or RFID) emit electromagnetic radiation, the function of the product may be impaired. High-frequency communication equipment can affect electrical medical products and impair their performance.

- ▶ Do not use the product in the vicinity of a magnetic resonance tomograph (MRT).
- ▶ Do not use the product next to or together with other devices. If such use is required, monitor the product and the other devices, and follow the relevant instructions for use in the event of malfunctions.
- ▶ Portable RF communications equipment including peripheral devices (e.g., antenna cables and external antennas) should be used no closer than 30 cm from the product, including cables specified by the manufacturer.
- ▶ Observe the information on electromagnetic compatibility; see chapter Electromagnetic compatibility [p. 26].
- ▶ In case of uncertainties, seek expert advice from KARL STORZ.
- ▶ Before use, a clinical/biomedical engineer or an EMC specialist should carry out an ad-hoc test of the electromagnetic radiation.

The use of accessories and cables other than those specified in the instruction manual may result in increased emissions or decreased immunity of the product. When using other accessories and cables, the operator is responsible for checking compliance with IEC 60601-1-2 for this particular product.

- ▶ To prevent increased electromagnetic emissions or reduced electromagnetic immunity of the product, only use accessories, transducers, and cables recommended or supplied by the manufacturer.

3.7 Patient leakage current

Patient leakage currents from products may add up if powered endoscopes and powered endotherapy devices are used simultaneously.

Excessively high leakage current levels may result in the patient becoming injured.

- ▶ Only use products of the same type together; e.g., CF.

3.8 National guidelines on airway management

Failure to observe the national guidelines on airway management may put the patient at risk.

- ▶ National guidelines on airway management must be observed in addition to the documentation accompanying the product.

3.9 Hot components

The high level of light intensity may cause the distal end, the light connections for the endoscope and fiber optic light cable, and adjacent components to heat up. This can cause burns to patients, users, or third parties.

The high level of light intensity may cause the distal end of the endoscope to heat up. This can cause burns to patients, users, or third parties.

- ▶ Set the light source output to a level that is just high enough to ensure optimal illumination of the operating area.
- ▶ Avoid contact with the distal end and light connections of the endoscope and fiber optic light cable.
- ▶ Avoid contact with the distal end of the endoscope.

3.10 High light intensity

The high level of light intensity produced by the light source may lead to permanent eye damage or blindness, and may cause tissue and items facing the light output to heat up.

- ▶ Do not look into the light output.
- ▶ Set the light source output to a level that is just high enough to ensure optimal illumination of the operating area.
- ▶ Make sure the light output is sufficiently far away from tissue and operating accessories.

4 Product description

4.1 Product overview



Video laryngoscope (8403AXC/BXC/DXC/GXC/KXC/MXC/HXP) without guidance for suction catheter

- 1 Socket for video connecting cable 2 Multifunction button



Video laryngoscope (8403EXC) without guidance for suction catheter



Video laryngoscope (8403AX/BX/HX) with guidance for suction catheter



Video connecting cable (8403X)

- 1 Laryngoscope connection socket 2 Monitor connection socket

4.2 Possible combinations

The C-MAC® Video Laryngoscopes from series 8403 can be operated with the following monitors:

- C-MAC® Monitor 8403ZX (software version: 704v300 or higher)
- C-MAC® Monitor 8404ZX
- C-MAC® Pocket Monitor 8403XD
- With other monitors via the camera control unit C-HUB® II (20290320) (software version: 721v105 or higher)

The C-MAC® Video Laryngoscopes from series 8403 can be combined with the following catheters:

- Blade 8403AX with 14 Fr. and 16 Fr. catheters
- Blade 8403HX and 8403BX with 16 Fr. and 18 Fr. catheters

4.3 Technical specifications

Video laryngoscope 8403AX/BX/HX/AXC/BXC/KXC/HXP

Designation	Value
Immersion protection	IPX8
Camera technology	CMOS
Resolution	640 x 480
Illumination:	LED, white, 1 W

Video laryngoscope 8403AX/BX/HX/AXC/BXC/KXC/HXP serial number 80.000 and higher













Designation	Value
Immersion protection	IPX8
Camera technology	CMOS
Resolution	300 x 400
Illumination	LED, white, 1 W

Video laryngoscope 8403DXC/EXC/GXC/MXC/NXC





Designation	Value
Immersion protection	IPX8
Camera technology	CMOS
Resolution	300 x 400
Illumination	LED, white, 1 W

4.4 Symbols employed

4.4.1 Symbols on the packaging

Symbol	Meaning
	Number of products in the product packaging
	Article no.
	CE conformity mark With this mark, the manufacturer declares the compliance of the products with the applicable regulation (EU) 2017/745. A code number after the CE mark indicates the responsible notified body.
	Consult instructions for use
	Manufacturer
	Serial number
	Keep dry
	Keep away from sunlight
	Temperature limit
	Date of manufacture
	In accordance with US federal law (21 CFR 801.109), this product may only be sold to or on prescription from a licensed physician.
	Not MR safe

4.4.2 Symbols on the product

Symbol	Meaning
	Follow instructions for use
	Type BF device
	Serial number
	CE conformity mark With this mark, the manufacturer declares the compliance of the products with the applicable regulation (EU) 2017/745. A code number after the CE mark indicates the responsible notified body.

4.5 Ambient conditions

Storage/transport conditions	
Temperature	-20°C ...+50°C
Relative humidity (non-condensing)	5–95%
Operating conditions	
Temperature	0°C ...40°C
Relative humidity (non-condensing)	30–70%
Air pressure	700–1,080 hPa

5 Preparation

5.1 Unpacking the product

1. Carefully remove the product and accessories from the packaging.
2. Check the delivery for missing items and evidence of shipping damage.
3. In the case of damage, hidden defects, and short deliveries, document their nature and extent and contact the manufacturer or supplier immediately.

5.2 Testing the product


5.2.1 Visual inspection

Missing components, surface changes, or other damage limit the life span of the product and may mean that the product can no longer be used for its intended use.

1. The product must be inspected for completeness, damage, and integrity before and after every use. Inspection tools such as magnifying glasses should be used as necessary.
2. Check for the following types of mechanical damage and changes:
 - Sharp corners
 - Burred edges
 - Rough surfaces
 - Protruding or bent parts



Defect pattern: Blade tip bent



Defect pattern: Blade roof bent


3. Check the product for rust and corrosion.

4. Check whether the coating is damaged.




Defect pattern: Poor general condition

5. Ensure that the components and connections are complete, intact, and positioned correctly.
6. Check whether there are residues and/or moisture at the interfaces.
7. Check the mobility of the components. Instruments that can be dismantled should be assembled for this purpose.
8. Use a magnifying glass to check that the telescope is complete, intact, and securely positioned.



Defect pattern: Upper glass missing



Defect pattern: Lower glass broken

5.2.2 Functional test

- ▶ Inspect the product for the following characteristics:
 - Fully functional control keys
 - Securely positioned connecting cable
 - Functional image transmission
 - Sufficient image quality
 - Light transmission

5.3 Connecting the C-MAC® Monitor

1. Insert the video connecting cable into the socket of the video laryngoscope. Pay attention to the orientation pin on the plug and socket.



2. Make sure that the cable is fully inserted.



3. Connect the other end of the video connecting cable to the socket on the C-MAC® Monitor.



4. Alternatively, connect the other end of the video connecting cable to the socket on the rear side.



- i** Another laryngoscope can be connected during use.

5.4 Connecting the C-MAC® Pocket Monitor

The Pocket Monitor can be connected to the product while folded up or down.

1. Connect the C-MAC® Pocket Monitor to the socket of the video laryngoscope. Pay attention to the orientation pin on the plug and socket.



2. Make sure that the lugs of the plug connection are opposite each other and the plug is fully inserted.
 - ⇒ As soon as the photo/video recording function is ready, the multifunction button lights up blue.
 - ⇒ The video laryngoscope is ready for operation.



i The monitor can be removed and connected to another laryngoscope during use.

5.5 Connecting C-Hub® II

1. Connect the video connecting cable to the socket on the product. Pay attention to the orientation pin on the plug and socket.



2. Make sure that the cable is fully inserted.



3. Connect the power supply unit to the C-HUB® II.
4. Connect the video connecting cable to the video input on the front of the C-HUB® II.
5. Connect the monitor to the S-Video output or HDMI output on the rear of the C-HUB® II.

6 Application

6.1 Using the product

⚠ WARNING

Wrong viewing mode! Risk of injury!

If the viewing mode settings are changed, this can alter the image display and image position.

- ▶ Before each use, check whether a live image or a saved image is being displayed and whether a correct image position has been selected.

Documentation can be carried out on the C-MAC® Monitor (8403ZX, 8404ZX) and C-MAC® Pocket Monitor (8403XD) using the multifunction button on the video laryngoscope.

1. Switch on the product in good time so that the optical system can warm up.
2. Switch on the monitor.
 - ⇒ As soon as the photo/video recording function is ready, the blue multifunction button lights up.



3. Press the multifunction button briefly to save images individually.



- ⇒ The multifunction button lights up green during the saving process.
4. Press the multifunction button for around 2 seconds to start video recording.
 - ⇒ During the video recording, the multifunction button lights up in pulsating green.
 5. Press the multifunction button briefly to stop video recording.

6.2 Inserting an oxygen or suction catheter

- ▶ Insert the catheter or the probe into the guide.



7 Disassembly

7.1 Disassembling the product

- ▶ Disconnect the C-MAC® Pocket Monitor or the connecting cable from the video laryngoscope.



8 Maintenance, servicing, repairs, and disposal

8.1 Repairs to the product

Repair work may only be performed by KARL STORZ or by a company authorized by KARL STORZ. The interventions described in this instruction manual are exempt from this rule.

- ▶ Please contact your local KARL STORZ subsidiary or authorized dealer (see the list of subsidiaries).

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, products must first be decontaminated. KARL STORZ reserves the right to send back contaminated products.

8.2 Disposing of the product

The product meets the requirements of the Directive on Waste Electrical and Electronic Equipment (WEEE).

Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this product.

1. The product must be disposed of in accordance with the applicable national laws and regulations at a suitable collection point for the reprocessing of electrical and electronic equipment.
2. Contact KARL STORZ SE & Co. KG, a KARL STORZ branch or an authorized dealer to find out the address of the collection point in your area.

9 Accessories and spare parts

9.1 Accessories

Item	Order no.
Forceps acc. to MAGILL, modified acc. to BOEDEKER, 25 cm	809125
Forceps acc. to MAGILL, for children, modified acc. to BOEDEKER, 20 cm	809120
Bag for intubation set C22, ULM model	8403YE
Protection cap for C-MAC® System Interface	8403YZ
C-MAC® Pocket Monitor connecting cable	8403XDP
C-MAC® Connecting Cable	8403X

10 Electromagnetic compatibility

10.1 Electromagnetic compatibility

The described product has been tested as a system with the following devices. Relevant electromagnetic compatibility (EMC) information can be found in the "Electromagnetic compatibility (EMC)" section of the instructions for use for the devices.

Device	Item number
C-MAC [®] Monitor	8403ZX
C-MAC [®] Monitor	8404ZX
C-MAC [®] Pocket Monitor	8403XD
C-HUB [®] II	20290320

The EMC warning statements, precautions, notes, and emission/immunity limits specified in the instructions for use for the devices also apply to the product described in these instructions for use.

11 Fault correction

11.1 Product malfunctions

Error description	Possible causes	Actions
Product failed	Power supply failure	▶ Check power supply
	Power plug and socket improperly connected	▶ Push the power plug firmly into the socket on the product
	System error	▶ Perform reset
No picture, TFT screen dark	Connecting cable not connected correctly	▶ Push the cable firmly into the socket on the product
	Defective camera electronics	▶ Send product to KARL STORZ for repair
	TFT screen defective	
Blurry picture, striping, streaking	Telescope is dirty	▶ Clean with cotton swab and alcohol solution or special cleaning paste
	Connecting cable contacts are dirty	▶ Clean the cable contacts
Color distortion	White balance not carried out correctly	▶ Repeat white balance adjustment
Color rendering alternates	Connecting cable broken	▶ Replace connecting cable
Image cannot be stored	No memory card inserted	▶ Insert memory card
	Memory card full	▶ Clean camera plug contacts
	Memory card is not recognized	▶ Insert the memory card while the monitor is switched off and then turn the monitor on again ▶ Perform reset ▶ Format the memory card on the PC with FAT
Video stream cannot be played back	No MPEG-4 Codec installed	▶ Install an MPEG-4 Codec
SD memory card cannot be written on	SD memory card is formatted incorrectly	▶ Format the SD memory card on the PC with FAT

12 Subsidiaries

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0, Fax: +49 7461 708-105
Email: info@karlstorz.com

KARL STORZ Endoskope Berlin GmbH

Scharnhorststr. 3, 10115 Berlin/Germany
Phone: +49 30 3069090, Fax: +49 30 30 19452

KARL STORZ Endoscopy Canada Ltd.

7171 Millcreek Drive, Mississauga, Ontario L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free (Canada only) Phone: 1-800-268-4880, Fax: 1-800-482-4198
(Canada only)
Email: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.

2151 East Grand Avenue, El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free (USA only) Phone: 800 421-0837, Fax: 800 321-1304 (USA only)
Email: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.

1 South Los Carneros Road, Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
Email: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.

815 N. W. 57th Avenue, Suite 480, Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
Email: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.

Edificio Atlantic, Oficina 3G, Calle D e/ 1ra y 3ra, 10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
Email: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.

Av. Ejercito Nacional No. 453 Piso 2, Colonia Granada, Alcaldia Miguel Hidalgo, C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
Email: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.

Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi, CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
Email: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.

Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez, Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
Email: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS

Stamveien1, 1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
Email: post@karlstorz.no

KARL STORZ Endoskop Sverige AB

Storsåtragränd 14, 127 39 Skårholmen, Sweden
Phone: +46 8 505 648 00
Email: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY

Taivaltie 5, 01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
Email: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG

Representation Office
Kęstučio st. 59 / Lenktoji st. 27, 08124 Vilnius, Lithuania
Phone: +370 5 272 0448, Mobile: +370 685 67 000
Email: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S

Skovlytoften 33, 2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
Email: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.

415 Perth Avenue, Slough, Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
Email: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.

Displayweg 2, 3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
Email: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.

Phone: +31 (0)33 4545890
Email: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.

12, rue Georges Guynemer, Quartier de l'Europe, 78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
Email: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH

Landstraßer Hauptstr. 148/1/G1, 1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
Email: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S. A.

Parque Empresarial San Fernando, Edificio Munich – Planta Baja, 28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
Email: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.

Via dell'Artigianato, 3, 37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
Email: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.

Capraška 6, 10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
Email: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.

Cesta v Gorice 34b, 1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
Email: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.

ul. Bojkowska 47, 44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
Email: info-pl@karlstorz.com

KARL STORZ Endoszkóp Magyarország Kft.

Toberek utca 2. fsz. 17/b, HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
Email: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl

Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4, 041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
Email: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.*

Patriarhou Grigoriou E' 34, 54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
Email: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**

Gedik Is Merkezi B Blok, Kat 5, D 38-39, Bagdat Cad. No: 162, Maltepe
Istanbul, Turkey

Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK

Derbenyevskaya nab. 7, building 4, 115114 Moscow, Russia

Phone: +7 495 983 02 40, Fax: +7 495 983 02 41

Email: info-ru@karlstorz.com

TOV LLC KARL STORZ Ukraine

Avenue Geroyiv Stalingrada Str. 2D, office 717 Kyiv, 04210/Ukraine

Phone: +38 095 000-895-0, +38-097-000-895-0, +38 073 000-895-0

Email: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG Representation Office

Sabit Orudschow 1184, apt. 23, 1025 Baku, Azerbaijan

Phone: +99 450 613 30 60

Email: info-az@karlstorz.com

KARL STORZ ENDOSKOPE – East Mediterranean and Gulf (Offshore) S.A.L.

Spark Tower 1st floor Charles Helou St., Horch Tabet – Sin El Fil, Beirut,
Lebanon

Phone: +961 1 501105, Fax: +961 1 501950

Email: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.

P.O. 6061, Roggebaai, 8012 Cape Town, South Africa

Phone: +27 21 417 2600, Fax: +27 21 421 5103

Email: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan

Saryarka, 6, BC "Arman", off. 910, 010000 Astana, Republic of Kazakhstan

Phone: +7 7172 552-549, 552-788, Fax: -444

Email: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf (branch)

Building West Side 7A – Unit 7WA – 3008, Dubai Airport Free Zone, P.O. Box
54983, Dubai - United Arab Emirates

Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282

Service Hotline: +971 (0)4 3415882

Email: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited

11th Floor, Dr. Gopal Das Bhawan, 28, Barakhamba Road, New Delhi
110001, India

Phone: +91 11 4374 3000, Fax: +91 11 4374 3010

Email: corporate@karlstorz.in

KARL STORZ SE & CO. KG

Interchange 21 Tower, Level 33, 399 Sukhumvit Road, North Klongtoey,
Wattana, 10110 Bangkok, Thailand

Phone: +84 28 3823 8000, Fax: +84 28 3823 8039

Email: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG

Resident Representative Office

14th Floor, MPlaza Saigon, 39 Le Duan, District 1, Ho Chi Minh City, Vietnam

Phone: +84 28 3823 8000, Fax: +84 28 3823 8039

Email: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.

Room 2503-05, 25F AXA Tower, Landmark East, No. 100 How Ming Street,
Kwun Tong, Kowloon, Hong Kong, People's Republic of China

Phone: +852 28 65 2411, Fax: +852 28 65 4114

Email: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch

Room 1805-1807, Building B, 18F Beijing IFC, No. 8, Jianguomenwai Street,
Chaoyang District, 100022, Beijing, People's Republic of China

Phone: +86 10 5638188, Fax: +86 10 5638199

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch

Room 701A Building 5 & Room 501 Building 7, No. 3000 Longdong Avenue,
Pilot Free Trade Zone, 201203, Shanghai, People's Republic of China

Phone: +86 21 60339888, Fax: +86 21 60339808

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch

Room 803-805, 8F Jin Jiang International Building, No. 1 West Linjiang
Road, Wuhou District, 6100414, Chengdu, People's Republic of China

Phone: +86 28 86587977, Fax: +86 28 86587975

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch

Room 2001-2005, 20F N-MEDIA International Center, No. 167 Youth Avenue,
Shenhe District, 110014, Shenyang, People's Republic of China

Phone: +86 24 23181118, Fax: +86 24 23181119

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch

Room 02B & 03 & 04A, 35F Teem Tower, No. 208 Tianhe Road, Tianhe
District, 510620, Guangzhou, People's Republic of China

Phone: +86 20 87321281, Fax: +86 20 87321286

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

Email: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

Email: infoasia@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Indonesia

Sinarmas MSIG Tower Level 37, Jl. Jend. Sudirman No. Kav. 21, Jakarta
Selatan

DKI Jakarta 12920

Email: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.

9F Hyowon-Building, 97, Jungdae-ro, Songpa-gu, 05719 Seoul, Korea

Phone: +82-70-4350-7474, Fax: +82-70-8277-3299

Email: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.

12F, No. 192, Sec. 2, Chung Hsin Rd., Sindian District, New Taipei City,
Taiwan

Phone: +886 933 014 160, Fax: +886 2 8672 6399

Email: info-tw@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Philippines

1901 Picadilly Star Bldg., 4th Avenue, BGC, Taguig City 1636, Philippines

Phone: +63 2 317 45 00, Fax: +63 2 317 45 11

Email: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.

Stage Bldg. 8F, 2-7-2 Fujimi, Chiyoda-ku, Tokyo 102-0071, Japan

Phone: +81 3 6380-8622, Fax: +81 3 6380-8633

Email: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd .

68 Waterloo Road, Macquarie Park NSW 2113, P O Box 50 Lane Cove NSW
1595, Australia

Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695

Toll free: 1800 996 562 (Australia only)

Email: info@karlstorz.au

www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE



USB825 • EN • V1.1 • 04-2021 • IFU • CE-MDR



KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-mail: info@karlstorz.com
www.karlstorz.com

